

HugeMed VLRM/VLDM

When time is critical, make your go-to device the HugeMed VLRM/VLDM.

Ideal for demanding OR, ED, and ICU settings, the VLRM/VLDM provides a consistently clear view of the airway, and enhanced maneuverability enabling fast intubations. The VLRM/VLDM system gives you a choice of 8 reusable blades and 5 disposable blades to cover a wide range of patients—from infant to morbidly obese.

Add real-time recording and an onboard tutorial, and you'll see why HugeMed VLRM/VLDM is the one you go to for advanced video laryngoscopy.

VLRM



Reusable Blades

VLM=ADVANCED AIRWAY VIEWS

- Digital color monitor with LCD—quality resolution
- Reusable video laryngoscope with digital-resolution camera
- Innovative blade sizing, angulation, and camera positioning
- Anti-fog feature, with a rapid heating profile to resist lens clouding and secretions

VLDM



Disposable Blades

VLM=ADVANCED USER FEATURES

- Integrated, real-time recording helps confirm tube placement and facilitate teaching
- Ergonomic design and rugged construction for multiple clinical applications
- Reusable blades in 8 sizes: MILLER 00, 0, 1, MAC 1,2,3,4,5
- Also available in 5 single-use sizes

Advanced Video Laryngoscopy:

Real-time recording. LCD clarity.



Applications

- First-use intubations, replacing direct laryngoscopy
- Normal or restricted oropharyngeal views
- Cormack–Lehane Grades I–IV laryngeal views
- Trauma airways—excellent when dealing with blood and secretions in the airway
- Video-guided tube exchange in the ICU



Specifications

Type	VLRM	VLDM
Display	8 Inch LCD monitor HD Display	
Resolution	1024*768RGB	
Data audiovisual output	Support	
Camera Resolution	2M Pixel	
View Angle	66°	
Illumination	≥800LUX	
Battery Life Time	> 2 hours	
Battery Capacity	3200mAh	
Voltage	3.7V	
Blade type	Reusable Blades	Disposable Blades
Material of Blade	SUS316L	PC
Shot and Record Function	Support	
Anti-fog	Anti-fog when it's on	
Blade Sizing	Miller 00/0/1; Mac 1/2/3/4/5	Mac 1/2/3/4/5
HDMI&USB Port	Include	
Connector for Screen and Handle	Include	
Bracket (with tray)	Include	
Attachments	Adapter, Power cord, Data connection line	
Blades Disinfection	Refer to the reusable laryngoscope instructions	Refer to the Disposable laryngoscope instructions